


## Pharmazie – Pharmaceuticals – Pharmacie



### Infusionsflaschen – Infusion bottles – Flacons perfusion


CA 2318	○	I	50,0	68,0	55	68,0	46,0	46,0	32,00
6328	○	II	50,0	68,0	55	68,0	46,0	46,0	32,00
6051-05	○	II	50,0	77,0	70	75,0	46,0	46,0	32,00
CA 2326	○	I	100,0	125,0	95	104,0	49,0	49,0	32,00
CA 2300	○	I	100,0	125,0	95	104,0	49,0	49,0	32,00
6198	○	II	100,0	128,0	95	104,0	49,0	49,0	32,00
6281	○	II	100,0	141,0	105	105,8	52,0	52,0	29,00
6281-01	○	II	100,0	141,0	105	105,8	52,0	52,0	29,00
6332	○	II	125,0	147,0	95	98,0	54,4	54,4	32,00
2321	○	I	150,0	207,0	139	102,5	66,0	66,0	32,00
6164	○	II	150,0	207,0	142	102,5	66,0	66,0	32,00
CA 2304	○	I	250,0	308,0	170	136,0	66,0	66,0	32,00

► Continued on page 8.2

## Pharmazie – Pharmaceuticals – Pharmacie



### Infusionsflaschen – Infusion bottles – Flacons perfusion



6199	○	II	250,0	308,0	170	136,0	66,0	66,0	32,00
6305	○	II	250,0	312,0	160	134,5	67,0	67,0	29,00
6305-01	○	II	250,0	312,0	160	134,5	67,0	67,0	29,00
6244-01	○	II	250,0	297,0	150	125,0	68,1	68,1	32,00
CA 2305	○	I	500,0	588,0	288	177,0	78,0	78,0	32,00
6184	○	II	500,0	588,0	225	177,0	78,0	78,0	32,00
6289	○	II	500,0	598,0	220	177,0	78,0	78,0	32,00
6282	○	II	500,0	582,0	260	177,0	78,6	78,6	29,00
6282-01	○	II	500,0	582,0	260	177,0	78,6	78,6	29,00
6245-01	○	II	500,0	570,0	220	147,0	86,0	86,0	32,00
6248	○	II	1000,0	1125,0	440	230,0	95,0	95,0	32,00
6248-01	○	II	1000,0	1125,0	440	230,0	95,0	95,0	32,00