



Tubular Glass Syringe systems



Bulk syringes Luercone/Luerlock

0.5 ml	47.6	6.85	4.65
1.0 ml long	54.0	8.15	6.35
1.0 ml standard	35.7	10.85	8.65
1.5 ml	43.2	10.85	8.65
2.25 ml	54.4	10.85	8.65
3.0 ml	72.2	10.85	8.65
5.0 ml	66.7	14.45	11.85



Material:

Type I glass (EP/USP/JP compliant), plunger rod PP, PS, Luerlock adapter PC

- Baked-on siliconisation as option (see page 6.1)
- Drug Master File Type III
- Suited to standard filling and packaging equipment
- Customised ceramic printing with heavy-metal-free inks
- Optional: special treatment with ammonium sulphate for greater surface resistance (see page 6.2)

Tip caps and rubber stoppers in various rubber formulations and designs available on request including the Gerresheimer TELC (Tamper-Evident Luerlock Closure).

Our technical support team would be happy to provide you with further information. Please contact info-buende@gerresheimer.com.