

Glass cartridges for diverse drug delivery systems

As an approved supplier of cartridges for various pharmaceutical applications our product range includes cartridges ranging from 1 to 10 ml.

To fulfil high and constantly increasing quality demands by our customers, all our cartridge lines are equipped with visual inspection systems for 100% dimensional control. In addition, our high quality lines make use of visual inspection systems to check for cosmetic defects and line scan cameras to control the glazed end.

Final packaging is performed under controlled environmental conditions.


Our product range includes clear and amber cartridges (glass type I) for the following applications:

- Cartridges for pen systems, auto- or needle-free injectors
- Insulin cartridges
- Dental cartridges

We offer ammonium sulphate treatment to improve the glass quality and processability (see page 6.2).



Tubular Glass Cartridges



Cartridges

61007	3.0	11.60	9.65	62.30
61001	3.0	11.60	9.70	62.30
62075	1.8	8.65	6.85	62.00
61006	1.5	8.65	6.85	56.95

Material:

Type I glass (EP/USP/JP compliant)

- Sizes from 1 to 10 ml, other sizes on request
- Packaging options: PP-box, tray, shrink wrap
- Suited to standard filling and packaging equipment
- Customised ceramic printing with heavy-metal-free inks
- Optional: special treatment with ammonium sulphate for greater surface resistance
- Cartridges (ISO or customised) to be used for pen systems, auto- or needle free injectors

Our technical support team would be happy to provide you with further information.
Please contact info-buende@gerresheimer.com.

