

### RTF® syringe accessories

We offer you a wide range of sterilised rubber stoppers suitable for all our RTF® syringe types. All rubber formulations are available as nested stoppers or in bags. Packaging in port bags for use with isolators or RABS is an additional option. For our 1 ml long syringes and 1 – 3 ml syringes, backstops are available on request.

Our range also includes plunger rods made from PP or PS for all syringe sizes from 0.5 ml – 5 ml. Customised versions are available on request.

Our technical support team will be happy to assist you in selecting the most appropriate combination. Please contact [info-buende@gerresheimer.com](mailto:info-buende@gerresheimer.com).

