



Tubular Glass Syringe systems



Readyject® syringe systems

0.5 ml	47.6	6.85	4.65
1.0 ml long	54.0	8.15	6.35
1.0 ml standard	35.7	10.85	8.65
2.25 ml	54.4	10.85	8.65
3.0 ml	72.2	10.85	8.65
5.0 ml	66.7	14.45	11.85

5/8 inch 25 G
1 inch 25 G
West® 4023/50 grey
Helvoet® FM 257/2



Material:

Type I glass (EP/USP/JP compliant), plunger rod PP, PS

- Readyject® as bulk version: glass barrel non-sterile; needle part gamma irradiated
- Readyject® as RTF®-version: glass barrel pre-sterilised; needle part gamma irradiated
- Tamper-evident rigid needle shield
- Baked-on siliconisation as option (see page 6.1)
- Drug Master File Type III
- Suited to standard filling and packaging equipment
- Customised ceramic printing with heavy-metal-free inks
- Optional: special treatment with ammonium sulphate for greater surface resistance (see page 6.2)

Our technical support team would be happy to provide you with further information.
Please contact info-buende@gerresheimer.com.