



Our next generation Autoinjector



Building bridges from
drugs to patients

- Single-use autoinjector
- Platform technology
- Patient-centric approach

gerresheimer
innovating for a better life

Single-use autoinjector for subcutaneous application

- Cartridge-based, pre-pressurized system
- Double needle system for shorter injection time
- No contact between the needle and the drug during storage
- Possibility of incorporating cartridges with baked-on silicone
- Terminal sterilization compatible with biologics



Remove from blister pack using pull tab

Platform technology enabling delivery of various drug formulations

- Adjustable spring force 30 to 70 N
- Adjustable subcutaneous thin wall needle 25, 27, and 29G
- Compatible with 1.5 and 3mL ISO cartridges



Patient – centric approach

- No release button, push-on-skin activation
- 2 simple user steps
- Partially transparent casing for visualisation of injection in real time



Push on skin activation

Simplicity of design

Cartridge based, pre-pressurized system

