



SensAIR™

Customizable on-body delivery device



Precise subcutaneous delivery
of large-volume biologics –
suitable for home use

The ready-to-use SensAIR™ on-body delivery device
can be customized to suit specific drugs and
patient needs:

- for a wide range of viscous liquids and flow rates
- automatic needle insertion and retraction
- small dimensions

The SensAIR™ on-body delivery device: safety and reliability through high precision

Gerresheimer has developed a novel pump technology for the delivery of biologics such as monoclonal antibodies. The wearable SensAIR™ on-body delivery device controls the subcutaneous administration of drugs of different viscosities. The device is designed to enable patients to self-administer up to 20 ml of biologics in a convenient way.

The SensAIR™ on-body delivery device design from Gerresheimer sets new standards in liquid drug delivery by:

Scalability

The design allows for a wide range of flow rates for biologics of different viscosities. The drug can be delivered from a 10 ml or a 20 ml cartridge.

Cost efficiency

The number of components and moving parts is kept to a minimum to help ensure cost efficiency.

Ready to use

The primary packaging is crucial for this ready-to-use prefilled and preloaded on-body delivery device. Gerresheimer provides both, the glass cartridge and the complete device.

Reliability

The SensAIR™ on-body delivery device is equipped with audible and visual signals. The injection can only start when the device is patched onto the skin and the needle is retracted automatically after use.

Accuracy

The innovative device architecture allows a controlled drug delivery of the cartridge content.

Connectivity

The optional Bluetooth® wireless connection to smartphones and computers enables easy connection of the device with mobile apps and other digital services.

Fully customized solutions to meet your requirements for the delivery of biologics



Device in development, not available for commercial use.

Are you interested in a reliable solution for the controlled subcutaneous administration of your biologic formulation?

Gerresheimer is your partner when it comes to state-of-the-art liquid drug delivery devices. Our platform of products enables cost-efficient customization to meet your needs.

Gerresheimer Olten
Phone: +41 62 209 71 00
info-olten@gerresheimer.com
gerresheimer.com

gerresheimer
innovating for a better life