



High-Quality Gx[®] Glass Ampoules for Parenteral Drugs



- ISO types B, C, D
- Color Break Ring
- Score Ring
- One Point Cut

gerresheimer.com

gerresheimer
innovating for a better life

Primary packaging tailored to your needs

Close partnership, technical support and comprehensive service

Gerresheimer is regarded as a technology leader and ranks among the leading companies in the world market. We offer a comprehensive range of ampoules from 1 to 30 ml made from borosilicate glass of hydrolytic class I

(Ph. Eur., USP, JP compliant). In addition to the ISO norms, diverse options enable us to match the type and properties of our ampoules to your special requirements and needs.

Colour Break Ring



Score Ring



One Point Cut



- More than 60 years experience in ampoules
- Ampoules made from flint or amber borosilicate glass of hydrolytic class I
- Sizes from 1 to 30 ml
- Straight-stem, funnel-type and closed ampoules (ISO types B, C and D) as well as customized special shapes
- Various break systems such as OPC (One Point Cut), CBR (Colour Break Ring) and SCR (Score Ring) with optimized opening forces
- Identification with up to three coloured code rings or customized ceramic printing with heavy-metal-free inks
- Complete inspection of the printing
- Special inner treatment with ammonium sulfate for greater surface resistance
- Inner siliconization for slide coating and optimum emptying
- Optimized designs for lyophilization
- State-of-the-art production technologies and camera inspection systems
- Excellent dimensional stability for high-speed filling and sealing operations
- Production and packaging under cGMP conditions
- DMF type III
- Production sites in Europe, Americas and Asia for risk mitigation

USA / Canada:

Gerresheimer Glass Inc. (Vineland)
Phone +1 856 692 3600
info-tg-us@gerresheimer.com

India: Gerresheimer

Pharmaceutical Packaging Mumbai Pvt. Ltd. (Kosamba Plant)
Phone +91 2646 270101 664
info-tg-in@gerresheimer.com

Europe:

Gerresheimer Wertheim
Phone +49 9342 9281-69
info-tg-euas@gerresheimer.com

Mexico /

Central and South America:

Gerresheimer Queretaro
Phone +52 442 441 1102
info-tg-us@gerresheimer.com

China: Gerresheimer

Shuangfeng Pharmaceutical Glass (Zhenjiang) Co. Ltd
Phone +86 511 888 980 00
info-tg-euas@gerresheimer.com

gerresheimer.com