





Primary packaging tailored to your needs Close partnership, technical support and comprehensive service

Gerresheimer is regarded as a technology leader and ranks among the leading companies in the world market. We offer a comprehensive range of vials from 0.6 to 50 ml in different designs, with or without blowback (European and American version) according to

industry standards or to your individual specifications. State-of-the-art production process and in-line camera inspection systems ensure high-performance in filling lines and lyophilization processes.

ISO vials



Screw thread vials



Double-chamber vials





- Expansion coefficients 33 and 51
- Sizes from 0.6 to 50 ml
- 13, 20, 32 mm finish, screw-thread finish, ISO and GPI standards
- Double entry screw thread capabilities
- Inner siliconization for optimum emptying
- Ammonium sulfate treatment to reduce alkali ion emission
- Optimized design for freeze drying processes
- Proprietary visual inspection systems ensure cosmetic and dimensional quality

- Standard, Pharma and PharmaPlus product line depending on your needs
- Production and packaging under cGMP conditions
- Ph. Eur., USP, JP, FEUM compliance
- DMF type III
- Production sites in Europe,
 Americas and Asia for risk mitigation
- Superior quality vials available:
 Gx® Armor, Gx® Elite and
 Gx® RTF (Ready To Fill)
- Customized vials formats to meet specific needs or new drug formats

Europe / Asia:

in vials

class I, II, III

Gerresheimer Boleslawiec

Phone +48 75 7350 432 info-tg-euas@gerresheimer.com

Product offering

- More than 60 years experience

- Vials made from flint or amber

borosilicate glass of hydrolytic

USA / Canada:

Gerresheimer Glass Inc. (Vineland)

Phone +1 856 692 3600 info-tg-us@gerresheimer.com

Mexico / Central and South America:

Gerresheimer Queretaro

Phone +52 442 441 1102 info-tg-us@gerresheimer.com

