



Perfectly Matched Moulded Glass Packaging

Close Partnership, Technical Support
and Comprehensive Service



Type I
Borosilicate Glass
Flint & Amber
2–600 ml

Type II
Sodalime Glass
with inner surface
treatment
Flint & Amber
5–50 ml

Type III
Sodalime Glass
Flint & Amber
Amber: 2–230 ml
Flint: 2–60 ml



Type I Offerings:
Injectible Range, ISO series,
DIN series, National, Victory,
Metro, Eros, Antibiotic, Infusion,
Boston, Wide Mouth Packer, Spray

Type III Offerings:
Injectible Range, ISO series,
Eros, Atlas, Galaxy, Poster Series,
Dropper Series, Syrup Range,
Essence Range, Wide Mouth Packer

Neck Profile:
Crimp, DIN, PP, Screw Neck

Body:
Graduated / Ungraduated

Applications:
Antibiotics, Veterinary, Oncology,
Surgeries, Essence, Ayurveda,
Food Supplements, Cosmetics,
Human Life Sciences

Production:
State of the art facility with
latest inspection technologies
ISO 9001, ISO 45001,
GMP Certified Facility



**Solution provider
with broad portfolio**

Registered Office:
906, Windfall Off. M. VasANJI Road
J.B. Nagar Andheri (East)
Mumbai 400 059 India

Manufacturing Sites:
Plot no. 38, Village Kunvarda
Taluka Mangrol Kosamba (RS)
District Surat - 394120

Telephone +91 22 2926 2360
support.india@gerresheimer.com

www.gerresheimer.com