



Tailor Made Solutions for Tubular Glass Packaging

Gerresheimer is a technology leader and ranks among leading companies in the world market offering comprehensive range



Vial Category Made of Type I Borosilicate in Amber & Flint Colors

High Quality HPLC Scientific Vials for Diagnostics

- Finishes: screw neck / crimp neck / press fit
- Add ons: Printed / Non Printed / Customized


High Quality Injectable Vials for Pharma & Parenteral use

- Sizes from 1-30 ml
- Finish: 13 mm / 20 mm;
- Screw neck / Crimp neck
- Add ons: Treated / Non Treated
- 100% Dimensional Gauging & Cosmetic Inspection



Ampoule Category: Made of Type I Borosilicate in Amber & Flint Colors

- Sizes from 1-30 ml
- Straight stem, Funnel type and Closed ampoules – ISO Form Type B, C, D
- Printed / Non Printed / Customised
- Identification upto 3 color code rings / customized ceramic printing



Color Break Ring Score Ring One Point Cut (OPC)



- Tubing options available in Type I Borosilicate
- State of the art production technologies and online camera inspection systems
- Excellent dimensional stability for high speed filling line and sealing operations
- Production and packaging under cGMP condition
- ISO 9001, ISO 45001 & GMP certified facility

Registered Office:
906, Windfall Off. M. VasANJI Road
J.B. Nagar Andheri (East)
Mumbai 400 059 India

Manufacturing Sites:
Plot no. 38, Village Kunvarda
Taluka Mangrol Kosamba (RS)
District Surat - 394120

Telephone +91 22 2926 2360
support.india@gerresheimer.com
www.gerresheimer.com