



Gerresheimer Primary Packaging Glass India

Close Partnership, Technical Support and Comprehensive Service
ISO 9001, ISO 45001 & GMP Certified Facilities

Moulded Glass India



Type I Borosilicate & Type III Sodalime Glass Vials & Bottles

Type I: Sizes from 2–600 ml in flint & amber

Offering for Injectable Range, ISO series, DIN series, National, Victory, Metro, Eros, Antibiotic, Infusion, Boston, Wide Mouth Packer

Neck Profile: Crimp, DIN, ISO, PP, Screw Neck
Body: Graduated / Ungraduated

Type III: Flint Sizes ranging from 2–60 ml & Type III Amber sizes ranging from 2–230 ml

Offering for Injectables, Eros, Atlas, Galaxy, Poster Series, Dropper Series, Syrup Range, Essence Range

Tubular Glass India



Vial Category Made of Type I Borosilicate in Amber & Flint Colors

High Quality HPLC Scientific Vials for Diagnostics

- Finishes: screw neck / crimp neck / press fit
- Add ons: Printed / Non Printed / Customised

High Quality Injectable Vials for Pharma & Parentals

- Sizes from 1–30 ml
- Add ons: Treated / Non Treated
- 100% Dimensional Gauging & Cosmetic Inspection

Ampoule Category Made of Type I Borosilicate in Amber & Flint Colors

- Sizes from 1–30 ml
- Straight stem, Funnel type and Closed ampoules – ISO Form Type B, C, D
- Add ons: Printed / Non Printed / Ceramic Printing / Customised
- Identification up to 3 color rings: Color Break Ring, Score Ring, One Point Cut (OPC)

Registered Office:
906, Windfall Off. M. VasANJI Road
J.B. Nagar Andheri (East)
Mumbai 400 059 India

Manufacturing Sites:
Plot no. 38, Village Kunvarda
Taluka Mangrol Kosamba (RS)
District Surat - 394120

Telephone +91 22 2926 2360
support.india@gerresheimer.com

www.gerresheimer.com