

GERRESHEIMER



High-Quality Gx[®] Glass Vials for Pharmaceuticals and Diagnostics



Primary packaging tailored to your needs

Close partnership, technical support and comprehensive service

Gerresheimer is regarded as a technology leader and ranks among the leading companies in the world market. We offer a comprehensive range of vials from 0.6 to 50 ml in different designs, with or without blowback (European and American version) according to industry standards or to your individual specifications. State-of-the-art production process and in-line camera inspection systems ensure high-performance in filling lines and lyophilization processes.

ISO vials



Screw thread vials



Double-chamber vials



PRODUCT OFFERING

- More than 60 years experience in vials
- Vials made from flint or amber borosilicate glass of hydrolytic class I, II, III
- Expansion coefficients 33 and 51
- Sizes from 0.6 to 50 ml
- 13, 20, 32 mm finish, screw-thread finish, ISO and GPI standards
- Double entry screw thread capabilities
- Inner siliconization for optimum emptying
- Ammonium sulfate treatment to reduce alkali ion emission
- Optimized design for freeze drying processes
- Proprietary visual inspection systems ensure cosmetic and dimensional quality
- Standard, Pharma and PharmaPlus product line depending on your needs
- Production and packaging under cGMP conditions
- Ph. Eur., USP, JP, FEUM compliance
- DMF type III
- Production sites in Europe, Americas and Asia for risk mitigation
- Superior quality vials available: Gx® Armor, Gx® Elite and Gx® RTF (Ready To Fill)
- Customized vials formats to meet specific needs or new drug formats



Europe / Asia: Gerresheimer Boleslawiec
Phone +48 75 7350 432 · info-tg-euas@gerresheimer.com

USA / Canada: Gerresheimer Glass Inc. (Vineland)
Phone +1 856 692 3600 · info-tg-us@gerresheimer.com

Mexico / Central and South America: Gerresheimer Queretaro
Phone +52 442 441 1102 · info-tg-us@gerresheimer.com

GERRESHEIMER

www.gerresheimer.com