

GERRESHEIMER



High-Quality Gx[®] Glass Vials
for Pharmaceuticals and
Diagnostics



Primary packaging tailored to your needs

Close partnership, technical support and comprehensive service

Gerresheimer is regarded as a technology leader and ranks among the leading primary packaging companies in the world market. With more than 60 years experience in tubular glass vials we offer a comprehensive range of flint and amber vials in different designs, with or without blowback (European and American version) according to industry standards or to your individual specifications. State-of-the-art production processes and in-line camera inspection systems ensure high-performance in filling lines and lyophilization processes.



PRODUCT OFFERING GERRESHEIMER INDIA

- Vials made from flint or amber borosilicate glass of hydrolytic class I
- Expansion coefficient 33 & 51
- Sizes from 2 to 30 ml
- Crimp neck, screw neck, printed and ammonium sulphate treated
- ISO 9001-2015 / ISO 15378-2015
- DMF Type III
- US, EP, JP compliant
- Proprietary visual inspection systems: Gx[®] RHOC for dimensional gauging and Gx[®] G3 for cosmetic inspection
- Other vial production sites in Europe and the Americas for risk mitigation

The brand new Kosamba Converting plant produces premium quality vials and ampoules.

It combines latest production technologies from the USA and Europe with a controlled environment throughout the plant. Commercial vial production to the highest international quality standards started at the end of 2016.

Gerresheimer Pharmaceutical Packaging Mumbai Private Limited

906, Windfall, Sahar Plaza Complex
Mathrudas Vasanji Road, J.B. Nagar, Andheri East
Maharashtra 400059 Andheri (East) Mumbai, India

Phone +919820192098
E-Mail s.sule@gerresheimer.com

GERRESHEIMER

www.gerresheimer.com